


SPECIFICATION



Mass Product

Contents

1. Description

1.1 General Description

1.2 Features

1.3 Application

1.4 Features Introduction

2. Product Specification

2.1 Size chart

2.1 Product Parameters

.application


Packing Criterion


. Handling Precautions

5.1 Product storage and maintenance

5.2 Handling Precautions








2. Product Specification

2.1 Size chart



size : 280*208*110mm

2.1 Product Parameters

Product Parameters		
1	Refond PN	RF-MTQAB7N03-B1
2	Battery	3.7V/2300mAh
3	Touch	RED+IR+UV) 20
		UV 20
		RED 20
4	LED	RED:650-665nm
		UV: 405-420nm
		IR: 850-860nm
5	LED Number of LED	203PCS
6	Charger	5V 2A
7		TYP-C
8	Power	8W
9		
10		330*250*140mm
11		680*520*300mm
12		8PCS

application



Mode 2: Blu light



Mode 3: red light



Packing Criterion

4.1 Package Diagram



Fig.3-2 Title




Fig.3-2 Title

. Handling Precautions

5.1 Product storage and maintenance

Do not expose the product to rain, humidity and high temperature.

Open flame sources, such as lit candles, should not be placed on the device

Careful supervision is required to prevent children from using this product.

If not used for a long time, charge the battery at least once every 3 weeks to maintain the battery life


Do not poke this product with sharp metal objects such as keys, hairpins, necklaces or nails.

5.2 Handling Precautions

Products for facial care, not recommended for body use.

Do not use this device for skin allergies, wounds, severe acne, burns, and bacterial skin infections.

Please consult a dermatologist before using this device if you are being treated for a skin disease or have a skin condition.



Date	Revisor	Version	Verifier	Remarks
2021-8-11		E/0		The first edition



Declare

This specification is written both in English and in Chinese and the latter is formal.